



PV Medical



4 IN 1 HIFU & LIPOHIFU SYSTEM

HIFU

Face Lift

Wrinkle Removal

Body Slimming



PV Medical




4 IN 1
HIFU & Lipohifu

4 in 1 face lift and body fat reduction machine

How lipohifu works?

Introduction

The Lipo hifu system uses high-intensity focused ultrasound deliver uniform energy at precise depth (13mm) to destroy targeted subcutaneous fat non-invasively, causing a rapid local temperature rise, destroying targeted adipose tissue just beneath the skin in the abdomen and flanks, predictable results after a single one-hour treatment. It can penetrate the skin layers and reach the targeted fat tissue without harming skin or surrounding tissues. Meanwhile the metabolized contents (triglyceride, fatty acid) of cells will excrete out of body automatically by blood circulation and lymphatic drainage without harming human body.



How HIFU Skin Tightening works?

Introduction

High Intensity Focused Ultrasound (HIFU) directly delivers heat energy to skin and subcutaneous tissue that can stimulate and renew the skin's collagen and thus consequently improving the texture and reducing sagging of the skin. It literally achieves the results of a facelift or a body lifts without any invasive surgery or injections, moreover, an added bonus of this procedure is that there is no downtime.


HIFU(High Intensity Focused Ultrasound) is delivered into the deep skin fractionally



Cause thermal coagulation to SMAS (Superficial Muscular Aponeurotic System), fascia, deep dermis and subcutaneous fibrous tissue.



By wound healing, skin tightening and rejuvenation happens gradually with time



12D HIFU




20,000 SHOTS
HIFU cartridges for
body and face




Single line

Multi line



600 SHOTS
Lipo HIFU cartridges
fat loss and body slimming



LIPOHIFU HANDLE



360° FACE HIFU CARTRIDGE

40,000 SHOTS

HIFU cartridges for
Face lift&compact



1.5mm




3.0mm



DEP INJECTION HANDLE

Vibration mode to increase blood circulation and comfort during operation



Non injection,
noninvasive
treatment



Effective introduction
of macromolecular
components



Instant introduction,
no harm to the skin



Vibration mode to
increase blood
circulation and comfort
during operation





Treatment range

Is a safe and effective means of non-invasive body sculpting.

- ✓ Parts all over the body
- ✓ Dissolve abdominal fat, reduce abdominal fat prominent
- ✓ Face, neck, arms and legs, buttocks and abdomen
- ✓ Postpartum and after liposuction reshape.
- ✓ Slimming, firming, improve cellulite and stretch marks

FACIAL AND BODY TREATMENT





“

Concerning the first experiences with FU4. Machine is working well. it is really enjoyable and impressive to customers to see the results on the screen.

”

Body parts ▼



Face ▼

