
ISO CE OEM ODM

Fractional RF Microneedle

Wrinkle Reduction

Skin Tightening

Acne Scars






Micro Needle RF Applicator

MICRO NEEDLE RF applicator for procedures requiring ablation and resurfacing of the skin.
Add vacuum suction combined probe for better skin contact, especially around eyes and neck areas.




NV Type

Non-Insulated Needle (frequency 2.2MHz)



IC Type

Insulated Needle (frequency 700KHz)




IS Type

Negative pressure non-insulated working head (output frequency 700KHz)



Fractional RF Appliicator

FRACTIONAL RF applicator for non-invasive wrinkle reduction suit for the patients fear of needle, no down time, recovery soon.








ICE Hammer


Provides cooling relief after microneedling RF treatments, reducing redness and swelling for a more comfortable recovery.

Face Treatment




- Non-surgical Face Lifting
- Wrinkle Reduction
- Skin Tightening
- Skin Rejuvenation
- Pore Reduction
- Acne Scars


Body Treatment





- Scars
- Stretch Marks














Principle introduction



Equipped with a negative pressure function, lift the skin



Non-insulated needles operate and continuously release RF energy.



Stimulate collagen regeneration



Before & After




Before




After



Before



After



Before



After



Before



After



Before



After



Before



After



Before



After



Before



After

Specification

Frequency	700KHz ± 20% / 2.2MHz ± 20%
Voltage	415Vpp / 110Vrms±10%
Power	45W (Max)
Vacuum Pressure	29mmHg (Max)
Vacuum Flow Rate	3.6L/min (Max)
Penetration Depth	7.0mm (Max)
Pressure	1.5kg (Max)
Power Supply	AC 110-220V 50/60Hz